



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁷ : A61J 1/05, 1/14	A1	(11) International Publication Number: WO 00/12043 (43) International Publication Date: 9 March 2000 (09.03.00)
<p>(21) International Application Number: PCT/SE99/01440</p> <p>(22) International Filing Date: 24 August 1999 (24.08.99)</p> <p>(30) Priority Data: 9802938-2 1 September 1998 (01.09.98) SE</p> <p>(71) Applicant (for all designated States except US): ASTRA AKTIEBOLAG [SE/SE]; S-151 85 Södertälje (SE).</p> <p>(72) Inventors; and (75) Inventors/Applicants (for US only): LUNDGREN, Anna [SE/SE]; Astra Hässle AB, S-431 83 Mölndal (SE). SUNDGREN, Mats [SE/SE]; Astra Hässle AB, S-431 83 Mölndal (SE).</p> <p>(74) Agent: ASTRA AKTIEBOLAG; Intellectual Property, Patents, S-151 85 Södertälje (SE).</p>	<p>(81) Designated States: AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CR, CU, CZ, DE, DK, DM, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZA, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SL, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).</p> <p>Published <i>With international search report. Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.</i></p>	
(54) Title: IMPROVED STABILITY FOR INJECTION SOLUTIONS		
<p>(57) Abstract</p> <p>A primary package containing a low molecular weight peptide-based thrombin inhibitors which package is sealed with a rubber stopper or plunger containing bromobutyl rubber.</p>		

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece	ML	Mali	TR	Turkey
BG	Bulgaria	HU	Hungary	MN	Mongolia	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MR	Mauritania	UA	Ukraine
BR	Brazil	IL	Israel	MW	Malawi	UG	Uganda
BY	Belarus	IS	Iceland	MX	Mexico	US	United States of America
CA	Canada	IT	Italy	NE	Niger	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NL	Netherlands	VN	Viet Nam
CG	Congo	KE	Kenya	NO	Norway	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NZ	New Zealand	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	PL	Poland		
CM	Cameroon	KR	Republic of Korea	PT	Portugal		
CN	China	KZ	Kazakhstan	RO	Romania		
CU	Cuba	LC	Saint Lucia	RU	Russian Federation		
CZ	Czech Republic	LI	Liechtenstein	SD	Sudan		
DE	Germany	LK	Sri Lanka	SE	Sweden		
DK	Denmark	LR	Liberia	SG	Singapore		
EE	Estonia						

IMPROVED STABILITY FOR INJECTION SOLUTIONS

Field of the invention

5 The present invention relates to solutions of low molecular weight thrombin inhibitors stored in primary packages containing rubber components, such as vials, bottles, cartridges and prefilled syringes. The invention also relates to the medical use of such stored thrombin inhibitor solutions.

10 Background of the invention

Solutions for parenteral use of pharmaceutically active substances are normally stored in primary packages such as, vials, bottles, cartridges or in prefilled syringes. The primary packages are sealed by a rubber stopper or plunger. A commonly used rubber material
15 contains chlorobutyl. Solutions of low molecular weight thrombin inhibitors stored in vials, bottles, cartridges and prefilled syringes sealed by a stopper or plunger containing chlorobutyl rubber exhibits increased degradation, leading to shortened time of storage.

Disclosure of the invention

20

It has now surprisingly been found that by using rubber material containing bromobutyl instead of chlorobutyl, the stability of the low molecular weight thrombin inhibitors in solution can be considerably improved.

25 The present invention provides a primary package, such as a vial, a bottle, a cartridge or a prefilled syringe containing a solution of a low molecular weight thrombin inhibitor for parenteral injection, sealed by a rubber stopper or plunger containing bromobutyl rubber instead of chlorobutyl rubber.

The present invention further provides a medical use of such thrombin inhibitor, or salts of such thrombin inhibitor, solutions kept in a primary package as mentioned above sealed by bromobutyl stoppers or plungers.

- 5 The present invention further provides an aqueous solution for parenteral administration comprising a low molecular weight peptide-based thrombin inhibitor or a salt thereof, having a pH in the range 3 to 8, preferably a pH about 5 and stored in a primary package, such as a vial, a bottle, a cartridge or a prefilled syringe, sealed by a rubber stopper or plunger containing bromobutyl.

10

- Thrombin inhibitors referred to in this application are low molecular weight peptide-based thrombin inhibitors. The term "low molecular weight peptide-based thrombin inhibitors" will be well understood by one skilled in the art to include thrombin inhibitors with one to four peptide linkages, and/or with a molecular weight below 1000, and includes those
15 described generically and, more preferably, specifically in the review paper by Claesson in Blood Coagul. Fibrin. (1994) 5, 411, as well as those disclosed in US Patent No. 4,346,078; International Patent Applications WO 97/23499, WO 97/02284, WO97/46577, WO 98/01422, WO 93/05069, WO93/11152, WO 95/23609, WO95/35309, WO 96/25426, WO 94/29336, WO WO 93/18060 and WO 95/01168; and European Patent Applications
20 623 596, 648 780, 468 231, 559 046, 641 779, 185 390, 526 877, 542 525, 195 212, 362 002, 364 344, 530 167, 293 881, 686 642, 669 317 and 601 459.

- Preferred low molecular weight peptide-based thrombin inhibitors include those known collectively as the "gatrans". Particular gatrans which may be mentioned include HOOC-
25 CH₂(R)Cha-Pic-Nag-H (known as inogatran; see International Patent Application WO 93/11152 and the list of abbreviations therein) and HOOC-CH₂-(R)Cgl-Aze-Pab-H (known as melagatran; see International Patent Application WO 94/29336 and the list of abbreviations therein).

30

The preferred low molecular weight peptide-based thrombin inhibitor to be kept in glass vials or syringes is selected from the group consisting of inogatran, (*Glycine, N*-[2-[2-[[[3-[(*aminoimino-methyl*)*amino*]propyl]*amino*]carbonyl]-1-piperidinyl]-1-(cyclohexylmethyl)-2-oxoethyl]-, [2*R*-[2*S*]]-), melagatran, (*Glycine, N*-[2-[2-[[[4(*aminoiminomethyl*)phenyl]-methyl]*amino*]carbonyl]-1-azetidiny]-1-cyclohexyl-2-oxoethyl]-, [2*R*-[2*S*]]-) and compound A, (*Glycine, N*-[1-cyclohexyl-2-[2-[[[4-[(*hydroxyimino*)aminomethyl]-phenyl]methyl]*amino*]carbonyl]-1-azetidiny]-2-oxoethyl]-, ethyl ester, [*S*-(*R*^{*}, *S*^{*})]-).

In one embodiment of the invention the thrombin inhibitor (preferably melagatran) solutions for parenteral injection is a water solution and are kept in primary packages such as vials, bottles, cartridges or prefilled syringes having a rubber stopper or plunger containing bromobutyl.

In another embodiment of the invention; the thrombin inhibitor for parenteral injection is in a water solution with an addition of hydroxy-propyl- β -cyclodextrin (HP β CD). The concentration of the thrombin inhibitor is in the range 0.001-100 mg/ml, preferably 2.5-20 mg/ml.

Working Example

Analytical technique

Liquid Chromatography (LC), for all analysis

The following equipment and parameters were used at the analysis of melagatran in solution.

Flowrate	1.0 ml/min
Wavelength	237 nm
Injection volume	20 µl
Analytical column	Waters Symmetry C8, 150 x 3.9 mm
5 Guard column	Waters Symmetry C8, 22 x 3.9 mm
Mobile phase	20 % (v/v) acetonitrile in phosphate buffer, pH 2.0 with 4.6 mM octanesulphonic acid.

EVALUATION

10

Results in tables are presented as total degradation of melagatran. This means that all by-products are included and presented as area% of melagatran.

Example 1.

15

This example shows a comparison of melagatran in HPβCD-solution in prefilled syringes (1.0 ml) having rubber plungers containing bromobutyl and chlorobutyl, respectively. The syringes were stored at 4, 25 and 50 °C for up to 6 months.

20 The melagatran solution was in direct contact with the different rubber materials.

MANUFACTURING OF SAMPLES

Melagatran, 2.5 mg/ml, in HPβCD water solution (40 % w/w), pH about 5

25 **Batch HF 839-2601**

Melagatran	442.1 mg
HPβCD	80.0 g
HCl, 1 M	qs
30 NaOH, 1 M	qs

water for injection to 200 g final weight (density 1.145 g/ml)

Melagatran was dissolved in water in a separate beaker and adjusted to pH 5.06. HP β CD powder was mixed with this solution together with water. The final solution was mixed
5 with a magnetic stirrer until the substance was completely dissolved and pH was finally adjusted to 5.02, and the solution was filtrated with a 0.22 μ m sterile filter.

Melagatran, 10 mg/ml, in HP β CD water solution (40 % w/w), pH about 5

10 **Batch HF 839-2602**

Melagatran	1.77 mg
HP β CD	80.0 g
HCl, 1 M	qs
15 NaOH, 1 M	qs
water for injection	to 200 g final weight (density 1.145 g/ml)

Melagatran was dissolved in water in a separate beaker and adjusted to pH 4.88. HP β CD powder was mixed with this solution together with water. The final solution was mixed
20 with a magnetic stirrer until the substance was completely dissolved and pH was finally adjusted to 5.0, and the solution was filtrated with a 0.22 μ m sterile filter.

FILLING OF SYRINGES (1.0 ml)

25 **Sample A1 (HF 839-2613) 10 mg/ml**

0.5 ml of HF 839-2602 was filled in 1 ml HYPACK[®] syringes from Becton Dickinson with a black plunger material (PH 701/50 from The West Company) containing chlorobutyl rubber.

Sample B1 (HF 839-2614) 10 mg/ml

0.5 ml of HF 839-2602 was filled in 1 ml HYPAK[®] syringes from Becton Dickinson with a grey plunger material (PH 4416/50 from The West Company) containing bromobutyl rubber.

Sample C1 (HF 839-2615) 2.5 mg/ml

0.5 ml of HF 839-2601 was filled in 1 ml HYPAK[®] syringes from Becton Dickinson with a grey plunger material (PH 4416/50 from The West Company) containing bromobutyl rubber.

Sample D1 (HF 839-2616) 10 mg/ml

0.5 ml of HF 839-2602 was filled in 1 ml HYPAK[®] syringes from Becton Dickinson with a black plunger material (PH 701/50 from The West Company) containing chlorobutyl rubber.

RESULTS OF STABILITY STUDIES**Sample A1 (HF 839-2613) 10 mg/ml - Chlorobutyl rubber**

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
0	5.2	-	1.2
1	5.2	4	1.0
1	5.3	50	7.4
3	5.1	4	1.2
3	5.1	25	4.5
3	5.2	50	14.9
6	5.1	4	1.2
6	5.1	25	3.7

Sample B1 (HF 839-2614) 10 mg/ml - Bromobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
0	5.2	-	1.1
1	5.2	4	1.0
1	5.2	50	6.4
3	5.1	4	1.2
3	5.1	25	2.4
3	5.2	50	12.8
6	5.1	4	1.1
6	5.1	25	3.1

Sample C1 (HF 839-2615) 2.5 mg/ml - Bromobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
0	5.3	-	1.2
1	5.4	4	1.1
1	5.3	50	7.2
3	5.3	4	1.3
3	5.3	25	3.9
3	5.2	50	14.2
6	5.2	4	1.2
6	5.2	25	5.7

Sample D1 (HF 839-2616) 10 mg/ml - Chlorobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
0	5.3	-	1.2
1	5.4	4	1.2
1	5.3	50	8.6
3	5.3	4	1.2
3	5.3	25	3.1
3	5.2	50	17.4
6	5.2	4	1.4
6	5.2	25	9.9

Conclusion

- 5 Rubber plungers containing chlorobutyl result in a more pronounced degradation compared to rubber plungers containing bromobutyl. This is true for high concentrations as well as low concentrations of melagatran in aqueous solutions.

10 The most pronounced difference was seen between plungers of chlorobutyl rubber and bromobutyl rubber when the dose of melagatran in aqueous solution was as low as 2.5 mg/ml.

Example 2.

- 15 This example is a comparison of melagatran in a water solution of HP β CD and melagatran in a water solution of NaCl. Both solutions are in direct contact with rubber plungers containing bromobutyl.

3 plungers of the quality FM 257 (from Helvoet Pharma N.V.) were placed in each 3 ml glass vial together with 1 ml solution of melagatran (NaCl water solution and HP β CD water solution, respectively). Reference samples, that is melagatran in NaCl water solution and in HP β CD water solution having no contact with plunger material. The reference
5 samples were treated in the same way as the other samples. The vials were stored at 50 °C for up to 3 months.

Compared to the study of Example 1 the ratio between solution exposed plunger surface and the quantity of melagatran solution is 16 times higher.

10

MANUFACTURING OF SAMPLES

Melagatran, 7.5 mg/ml, in HP β CD water solution (40 % w/w), pH about 5.

Batch HF 839-2679

15	Melagatran	928.8 mg
	HP β CD	55.0 g
	HCl, 1 M	qs
	NaOH, 1 M	qs
	water for injection	137.4 g (density 1.145 g/ml)

20

Melagatran and HP β CD were dissolved in water and adjusted to pH 4.96. The final solution was diluted with water to final weight and sterile filtrated with 0.45 μ m filter.

25 **Melagatran, 7.5 mg/ml, in NaCl water solution, pH about 5.**

Batch HF 839-2680

	Melagatran	1315.5g
	NaCl	1.441 g
30	HCl, 1 M	qs

NaOH, 1 M	qs
water for injection	to 170 (density 1.0 g/ml)

Melagatran and NaCl were dissolved in water and adjusted to pH 5.03. The final solution
5 was diluted with water to final weight and sterile filtrated with 0.22 µm filter.

FILLING OF VIALS

Sample A2 (HF 839-2682) 7.5 mg/ml in NaCl

10 1.0 ml of HF 839-2680 was filled in 3 ml vials together with 3 black unsiliconized plungers
(FM 257 from Helvoet Pharma N.V.) containing bromobutyl rubber.

Sample B2 (HF 839-2683) 7.5 mg/ml in NaCl

1.0 ml of HF 839-2680 was filled in 3 ml vials together with 3 black siliconized plungers
15 (FM 257 from Helvoet Pharma N.V.) containing bromobutyl rubber.

Sample C2 (HF 839-2684) 7.5 mg/ml in NaCl

1.0 ml of HF 839-2680 was filled in 3 ml vials together with 3 grey siliconized plungers (FM 257 from Helvoet Pharma N.V.) containing bromobutyl rubber.

Sample D2 (HF 839-2688) 7.5 mg/ml in NaCl

1.0 ml of HF 839-2680 was filled in 3 ml vials (Reference).

Sample E2 (HF 839-2689) 7.5 mg/ml in HP β CD

25 1.0 ml of HF 839-2679 was filled in 3 ml vials together with 3 black unsiliconized plungers (FM 257 from Helvoet Pharma N.V.) containing bromobutyl rubber.

Sample F2 (HF 839-2690) 7.5 mg/ml in HP β CD

1.0 ml of HF 839-2679 was filled in 3 ml vials together with 3 black siliconized plungers (FM 257 from Helvoet Pharma N.V.) containing bromobutyl rubber.

Sample G2 (HF 839-2691) 7.5 mg/ml in HP β CD

1.0 ml of HF 839-2679 was filled in 3 ml vials together with 3 grey siliconized plungers (FM 257 from Helvoet Pharma N.V.) containing bromobutyl rubber.

Sample H2 (HF 839-2695) 7.5 mg/ml in HP β CD

1.0 ml of HF 839-2679 was filled in 3 ml vials (Reference).

RESULTS OF STABILITY STUDIES**Sample A2 (HF 839-2682) 7.5 mg/ml in NaCl - Bromobutyl rubber**

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
1	5.9	50	4.2
3	6.0	50	9.3

Sample B2 (HF 839-2683) 7.5 mg/ml in NaCl - Bromobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
1	5.8	50	4.0
3	6.0	50	8.7

Sample C2 (HF 839-2684) 7.5 mg/ml in NaCl - Bromobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
1	5.8	50	3.7
3	5.8	50	7.9

Sample D2 (HF 839-2688) 7.5 mg/ml in NaCl - Reference

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
1	5.2	4	1.4
3	5.3	4	1.4
1	5.4	50	3.4
3	5.6	50	6.8

Sample E2 (HF 839-2689) 7.5 mg/ml in HP β CD - Bromobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % melagatran)
1	5.5	50	5.5
3	5.6	50	11.3

5

Sample F2 (HF 839-2690) 7.5 mg/ml in HP β CD - Bromobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
1	5.4	50	5.4
3	5.5	50	11.3

Sample G2 (HF 839-2691) 7.5 mg/ml in HP β CD - Bromobutyl rubber

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
1	5.4	50	5.4
3	5.5	50	10.3

Sample H2 (HF 839-2695) 7.5 mg/ml in HP β CD - Reference

Storage time (months)	pH	Temperature (°C)	Total degradation (area % of melagatran)
1	5.2	4	1.5
3	5.3	4	1.7
1	5.3	50	5.7
3	5.4	50	10.7

Conclusion

5

Melagatran in a water solution of NaCl exhibits a somewhat lower degradation compared to melagatran in a water solution of HP β CD. This is true both for solutions in contact with plunger material (FM 257 bromobutyl) 8%* compared to 11%*, and solutions in absence of plunger material (reference) 7%* compared to 11%*.

10

*; is total degradation in area% of melagatran

Example 3.

15

This example shows a comparison of different kinds of stopper and plunger materials containing either bromobutyl rubber or chlorobutyl rubber in contact with a melagatran solution (NaCl, pH 5). Melagatran solution was filled in glass vials (3 ml) together with stoppers and plungers of different brands. 5 different rubber materials were used in the study. There were 3 different bromobutyl and 2 different chlorobutyl rubbers. As reference,

20

NaCl water solution of melagatran was stored without any contact with stopper or plunger material.

The ratio between exposed plunger or stopper surface and melagatran in water solution is higher than in Example 1. A calculation has been made of exposed area of each tested plunger or stopper material. In the study the area ratio is 10-15 times higher compared to the area represented in Example 1. The vials were studied up to 19 days at a temperature of 50°C.

MANUFACTURING OF SAMPLES

Melagatran, 5 mg/ml, in isotonic NaCl solution, pH about 5.

Batch HF 839-2719

Melagatran	10.0 mg
NaCl	17.6 g
HCl, 1 M	qs
NaOH, 1 M	qs
water for injection	To 2000 g final weight (density 1.0 g/ml)

Melagatran and NaCl were dissolved in water and pH adjusted to 4.95 The solution was diluted to final weight with water.

FILLING OF VIALS

The total contact surface between the rubber material and the solution was enhanced in different ways and different extent. One way was by putting pieces of vial stopper material into each vial. For sample A3, the stopper material was divided into eight equal parts, and two parts in each vial (total of 2/8). Another way to enhance the contact surface was to put 2-3 plungers in each vial. For sample E3, three plungers were put in each vial. In samples A3 to F3, the contact surface was increased of 10-15 times compared to the normal contact surface between plunger and solution in a 1 ml syringe (used in Example 1).

Sample A3 (HF 839-2727) 5 mg/ml in NaCl

1.5 ml of HF 839-2719 was filled in a 3 ml vial together with two 1/8 parts of a 10 ml vial stopper (FM 50 from Helvoet Pharma N.V.) containing chlorobutyl rubber.

5

Sample B3 (HF 839-2728) 5 mg/ml in NaCl

1.5 ml of HF 839-2719 was filled in 3 ml vial together with 2 grey plungers (PH 4023/50 from The West Company) containing bromobutyl rubber.

10

Sample C3 (HF 839-2729) 5 mg/ml in NaCl

1.5 ml of HF 839-2719 was filled in 3 ml vial together with 2 black plungers (PH 701/50 from The West Company) containing chlorobutyl rubber.

Sample D3 (HF 839-2730) 5 mg/ml in NaCl

15 1.5 ml of HF 839-2719 was filled in 3 ml vial together with 2 grey plungers (W 4416/50 from The West Company) containing bromobutyl rubber.

Sample E3 (HF 839-2731) 5 mg/ml in NaCl

20 1.5 ml of HF 839-2719 was filled in 3 ml vial together with 3 black plungers (FM 257 from Helvoet Pharma N.V.) containing bromobutyl rubber.

Sample F3 (HF 839-2732) 5 mg/ml in NaCl

1.5 ml of HF 839-2719 was filled in 3 ml vial (Reference).

RESULTS OF STABILITY STUDIES

Sample A3 (HF 839-2727) 5 mg/ml in NaCl - Chlorobutyl rubber

Storage time (days)	pH	Temperature (°C)	Total degradation (area % of melagatran)
11	~5.0	50	8.0
19	~5.0	50	11.8

5 **Sample B3 (HF 839-2728) 5 mg/ml in NaCl - Bromobutyl rubber**

Storage time (days)	pH	Temperature (°C)	Total degradation (area % of melagatran)
11	~5.0	50	0.9
19	~5.0	50	1.4

Sample C3 (HF 839-2729) 5 mg/ml in NaCl - Chlorobutyl rubber

Storage time (days)	pH	Temperature (°C)	Total degradation (area % of melagatran)
11	~5.0	50	1.5
19	~5.0	50	2.4

Sample D3 (HF 839-2730) 5 mg/ml in NaCl - Bromobutyl rubber

Storage time (days)	pH	Temperature (°C)	Total degradation (area % of melagatran)
11	~5.0	50	1.3
19	~5.0	50	1.6

Sample E3 (HF 839-2731) 5 mg/ml in NaCl - Bromobutyl rubber

Storage time (days)	pH	Temperature (°C)	Total degradation (area % of melagatran)
11	~5.0	50	1.2
19	~5.0	50	1.4

Sample F3 (HF 839-2732) 5 mg/ml in NaCl - Reference

Storage time (days)	pH	Temperature (°C)	Total degradation (area % of melagatran)
11	~5.0	50	0.6
19	~5.0	50	1.0

Conclusion

All three bromobutyl rubber materials demonstrate lower melagatran degradation compared to the two chlorobutyl rubber materials.

Summary conclusion

It is shown in Example 1 that, for water solutions containing melagatran stored in HYPAK[®] syringes (from Becton Dickinson), improved stability is demonstrated using plungers containing bromobutyl rubber compared to the corresponding plungers containing chlorobutyl rubber.

It is shown in Example 2 that, for water solutions of melagatran stored in glass vials, improved stability is demonstrated using a NaCl water solution compared to a HP β CD water solution. This is true for melagatran in solution with and without contact of plungers containing bromobutyl rubber.

It is shown in Example 3 that for melagatran in a NaCl water solution, improved stability is demonstrated using rubber materials containing bromobutyl compared to rubber materials containing chlorobutyl.

CLAIMS

1. A primary package containing an aqueous solution for parenteral administration
5 comprising a low molecular weight peptide-based thrombin inhibitor or a salt thereof,
having a pH in the range 3 to 8, the primary package being sealed with a rubber stopper or
plunger containing bromobutyl rubber.
2. A primary package according to claim 1, wherein the primary package is a vial.
- 10 3. A primary package according to claim 1, wherein the primary package is a bottle.
4. A primary package according to claim 1, wherein the primary package is a cartridge.
- 15 5. A primary package according to claim 1, wherein the primary package is a prefilled
syringe.
6. A primary package according to any of the preceding claims, wherein the solution is a
NaCl solution.
- 20 7. A primary package according to any of the preceding claims, wherein the solution also
comprises hydroxy-propyl- β -cyclodextrin.
8. A primary package according to any of the preceding claims, wherein the concentration
25 of the thrombin inhibitor in the solution is in the range 0.001-100 mg/ml, preferably 2.5-20
mg/ml.
9. A primary package according to any of the preceding claims, wherein the pH of the
solution is in the range 3-8.

10. A primary package according to claim 9, wherein the pH of the solution is about 5.
11. A primary package according to any of the preceding claims, wherein the thrombin inhibitor is melagatran.
- 5 12. A primary package according to any of the preceding claims, wherein the thrombin inhibitor in the solution is inogatran.
13. A primary package according to any of the preceding claims, wherein the thrombin
10 inhibitor in the solution is compound A.
14. A primary package according to any of the preceding claims, wherein the bromobutyl rubber material consists of, or correspond to, the quality PH 4023/53.
- 15 15. A primary package according to any of claims 1-13, wherein the bromobutyl rubber material consists of, or correspond to, the quality W 4416/50.
16. A primary package according to any of claims 1-13, wherein the bromobutyl rubber material consists of, or correspond to, the quality FM 257.
- 20 17. Use of a rubber stopper or plunger containing bromobutyl rubber for sealing a primary package, such as a vial, a bottle, a cartridge or a prefilled syringe, containing a low molecular weight peptide-based thrombin inhibitor in an aqueous solution.
- 25 18. A process for the manufacture of a primary package according to claim 1 comprising the steps of dissolving a low molecular weight peptide-based thrombin inhibitor in an aqueous solution, adjusting the pH of the solution to be in the range 3 to 8, optionally adding a cyclodextrin substance, sterile filtering the solution and filling it on a primary package which is then sealed with a rubber stopper or plunger containing
30 bromobutyl rubber.

INTERNATIONAL SEARCH REPORT

International application No.

PCT/SE 99/01440

A. CLASSIFICATION OF SUBJECT MATTER

IPC7: A61J 1/05, A61J 1/14

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC7: A61J

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

SE,DK,FI,NO classes as above

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	WO 9633216 A1 (PHARMACIA AB), 24 October 1996 (24.10.96), page 6, line 17 - line 19 --	1-18
X	EP 0390244 A1 (DUPHAR INTERNATIONAL RESEARCH B.V.), 3 October 1990 (03.10.90), page 2, line 48 - line 51; page 3, line 14 - line 15, page 4, Table A --	1-18
A	WO 9739770 A1 (ASTRA AKTIEBOLAG), 30 October 1997 (30.10.97) -- -----	1-18

☐ Further documents are listed in the continuation of Box C.☒ See patent family annex.

* Special categories of cited documents:

- "A" document defining the general state of the art which is not considered to be of particular relevance
- "E" earlier document but published on or after the international filing date
- "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- "O" document referring to an oral disclosure, use, exhibition or other means
- "P" document published prior to the international filing date but later than the priority date claimed

"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

"X" document of particular relevance: the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone

"Y" document of particular relevance: the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art

"&" document member of the same patent family

Date of the actual completion of the international search

Date of mailing of the international search report

27 December 1999

Name and mailing address of the ISA/
Swedish Patent Office
Box 5055, S-102 42 STOCKHOLM
Facsimile No. +46 8 666 02 86

Authorized officer

Nebil Gecer/Els
Telephone No. +46 8 782 25 00

INTERNATIONAL SEARCH REPORT
Information on patent family members

02/12/99

International application No.

PCT/SE 99/01440

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 9633216 A1	24/10/96	AU 5412996 A	07/11/96
		EP 0871959 A	21/10/98
		IL 117907 D	00/00/00
		JP 10510925 T	20/10/98
		SE 9501472 D	00/00/00
		US 5862196 A	19/01/99
EP 0390244 A1	03/10/90	SE 0390244 T3	
		AT 84205 T	15/01/93
		CA 2012919 A	28/09/90
		DK 390244 T	08/02/93
		IL 93881 A	07/10/94
		JP 2283377 A	20/11/90
		NL 8900747 A	16/10/90
		US 5383864 A	24/01/95
WO 9739770 A1	30/10/97	AU 2719597 A	12/11/97
		AU 7716496 A	19/06/97
		EP 0864173 A	16/09/98
		SE 9601556 D	00/00/00